

Recombinant Human Bone Morphogenetic Protein-7 (BMP-7)

Certificate of Analysis and Data Sheet

➤ Source: E.Coli	➤ Catalog No. CTK-333
----------------------------	---------------------------------

➤ **Background :**

BMP-7 (osteogenic protein 1) belongs to the Bone-growth regulatory factors that are members of the transforming growth factor-beta (TGF-beta) superfamily of proteins. They are synthesized as large precursor molecules which are cleaved by proteolytic enzymes. The active form can consist of a dimer of two identical proteins or a heterodimer of two related bone morphogenetic proteins.

➤ **Description :**

Recombinant Human BMP-7 produced in E.Coli is a **monomeric**, non-glycosylated, Polypeptide chain containing 139 amino acids and having a molecular mass of 15679.97 Dalton.

BMP-7 is purified by proprietary chromatographic techniques.

➤ **Physical Appearance:**

Sterile Filtered White Lyophilized (freeze-dried) powder.

➤ **Formulation:**

BMP-7 was lyophilized from a concentrated (1mg/ml) sterile solution containing 10mM sodium citrate pH=3.5.

➤ **Solubility:**

It is recommended to reconstitute the lyophilized BMP-7 in sterile 18MΩ-cm H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

➤ **Stability:**

Lyophilized BMP-7 although stable at room temperature for 3 weeks, should be stored desiccated below -18 C. Upon reconstitution BMP-7 should be stored at 4 C between 2-7 days and for future use below -18 C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).

Please avoid freeze-thaw cycles.

➤ **Purity:**

Greater than 95.0% as determined by:

(a) Analysis by RP-HPLC.

(b) Anion-exchange FPLC.

(c) Analysis by reducing and non-reducing SDS-PAGE Silver Stained gel.

➤ **Amino acid sequence:**

The sequence of the first five N-terminal amino acids was determined and was found to be Ser-Thr-Gly-Ser-Lys.

➤ **Dimers and aggregates:**

Less than 1% as determined by silver-stained SDS-PAGE gel analysis.

➤ **Applications:**

1. Molecular standard (Western, ELISA) in studying secreted BMP-7
2. Preparing antibodies for BMP-7 monomer.
3. Molecule standard in detecting secreted BMP-7 in reduced SDS-PAGE.

➤ **Endotoxin:**

Less than 0.1 ng/μg (IEU/μg) of Bone Morphogenetic protein-7 .

➤ **Protein content:**

Protein quantitation was carried out by two independent methods:

1. UV spectroscopy at 280 nm
2. Analysis by RP-HPLC, using a standard solution of BMP-7 as a Reference Standard.

➤ **Usage:**

This material is offered for research, laboratory or further evaluation purposes.